

PRODUCT DESCRIPTION - PD 278509-1.0EN

Axtra® PHY GOLD 30 T

Phytase feed enzyme

Application

Axtra® PHY GOLD 30 T is a phytase feed enzyme system specifically developed to increase the digestibility of phytin-bound phosphorus, calcium, energy and amino acids in pig and poultry diets.

Usage levels

Use at a rate of 0.008-0.13 kg/tonne (0.0008-0.013%) of finished feed, included either directly or via a premix.

Processing stability

For optimum bioefficacy, do not exceed conditioning and pelleting temperatures of 95°C (203°F).

Composition

Dried *Trichoderma reesei* fermentation product
Sodium sulphate
Poly-vinyl alcohol (PVA)
Hydrated magnesium silicate

Guaranteed minimum activity:

6-phytase (EC 3.1.3.26) 30000 FTU/g

Analysis methods detailing assay conditions are available upon request.

Properties

Axtra® PHY GOLD 30 T is supplied as an off-white to light tan, fine granular product as follows:

Typical particle size based on sieve fineness:

Micrometers	Approx %
< 250	5
250 to 500	90
> 500	5

Particle size distribution studies by Method DIN 53734 indicate absence of particles less than 10 micrometers.

Axtra® PHY GOLD 30 T is sieved to 1.2 mm prior to packaging.

Bulk density	1000 - 1500 kg/m ³
Moisture content	Maximum 2 %

Physical/chemical specifications

Axtra® PHY GOLD 30 T is produced in accordance with FAO/WHO JECFA and FCC recommendations as follows:

Lead	< 5 ppm
------	---------

Microbiological specifications

Axtra® PHY GOLD 30 T is produced in accordance with FAO/WHO JECFA and FCC recommendations as follows:

Aflatoxins	Negative by test
Antibiotic activity	Negative by test
Salmonella/25g	Negative by test
E.coli/25g	Negative by test
Coliforms	< 30 CFUs/g

There are no producing micro-organisms in the final product. None of the ingredients is of animal origin.

Storage

Store in dry conditions. The constituent enzyme activities remain stable for at least 12 months when stored in original packaging at <25°C (<77°F) and 6 months when included in a vitamin/mineral premix and stored at <25°C (<77°F).

PRODUCT DESCRIPTION - PD 278509-1.0EN

Axtra® PHY GOLD 30 T

Phytase feed enzyme

Packaging

Axtra® PHY GOLD 30 T is available in 25kg multi-wall polyethylene lined paper bags and in 1000kg bulk (tote) bags.

Safety and handling

Provide for good ventilation and avoid dust formation. It is recommended to use protective glasses, respiratory mask and gloves during handling. In case of accidental contact with skin or eyes, the only action needed is copious flushing with water. See the Material Safety Data Sheet for further information.

Additional information

For further information please contact:
DuPont Animal Nutrition

EMEA: Genencor International, B.V.
Tel +31 71 5686 168

NAM: Danisco US Inc.
Tel +1 800-847-5311

LATAM: Danisco Argentina S.A.
Tel +54 357 6429000

ASPAC: Danisco Singapore Pte Ltd.
Tel +65-6586-3688

China: Genencor (China) Bio-Product Co., Ltd.
Tel +86 510 8370 9751

<http://animalnutrition.dupont.com/>
info.animalnutrition@dupont.com